Gelatin Sealed ePTFE Vascular Prostheses

Instructions for Use

English	Instructions For Use	5
Français	Instructions D'Utilisation	10
Deutsch	Gebrauchsanweisung	16
Nederlands	Gebruiksaanwijzing	22
Italiano	Istruzioni Per L'Uso	28
Español	Instrucciones De Empleo	34
Português	Instruções De Utilização	40
Svensk	Bruksanvisning	46
日本語	取扱説明書	51
Dansk	Brugsvejledning	56
Norsk	Brukerveiledning	61
ΕΛΛΗΝΙΚΑ	Οδηγίες Χρήσης	66
Česky	Návod k použití	72
Magyar	Használati útmutató	77
Lietuvių Kalba	Naudojimo instrukcijos	82
Slovensky	Návod na použitie	87
Polska wersja	Instrukcja dla użytkownika	92
中文	使用说明	98
Русский	Инструкции по применению	102
Türkçe	Kullanım Talimatları	108

English Instructions for Use

Caution: Federal (USA) law restricts this device to sale, distribution and use by or on the order of a physician.

Description

Vascutek Ltd gelatin sealed ePTFE grafts are manufactured in both supported and unsupported ePTFE (expanded polytetrafluoroethylene) and sealed with absorbable gelatin on the external surface only. The gelatin is the same as that used on Vascutek Ltd sealed polyester grafts¹ and has been shown to hydrolyze within approximately 14 days². The external support, where applicable, provides kink resistance. The end supported and fully supported grafts are designed to be peeled to facilitate the fashioning of the anastomoses. The central spiral support is **NOT** designed to be removed. Wrapped grafts have an outer layer of ePTFE, which is **NOT** designed to be removed. All grafts are supplied sterile and will remain so for the stated shelf life, unless the packaging is opened or damaged.

Origin Of Gelatin

Vascutek Ltd uses gelatin manufactured from animals native to and exclusively raised in Australia. Australia is one of only a few countries recognised as free from TSE infected animals, including BSE and scrapie. The EU Scientific Steering Committee has conducted Geographical BSE Risk Assessment (GBR) and concluded that Australia has the most favourable level 1 rating in relation to BSE risk.

Indications

Vascutek Ltd gelatin sealed ePTFE grafts are used for systemic vascular repair (excluding Central Circulatory System as defined by Annex IX, 1.7 of the Medical Device Directive 93/42/EEC) except for the application of the Blalock-Taussig procedure³, examples of which are indicated in the table below. In Canada, Vascutek Ltd. ePTFE grafts are intended for the creation of subcutaneous arteriovenous conduits for blood access, bypass or reconstruction of occluded or diseased arterial blood vessels.

Configuration	Graft Characteristic		Typical Applications
	Supported	Unsupported	
Straight Standard Wall	High Strength & Compression Resistant	-----	Systemic vascular repair, primarily for axillo-femoral/ bi-femoral bypass and femoral popliteal reconstruction.
Straight Standard Wall	-----	High Strength	Systemic vascular repair, primarily for femoral popliteal reconstruction and the creation of arteriovenous conduits for blood access.
Straight - Central Support Standard Wall	High Strength & Compression Resistant	-----	Arteriovenous conduit for blood access.
Straight - End Support Standard Wall	High Strength & Compression Resistant	-----	Arteriovenous conduit for blood access.
Straight Thin Wall	Thin, Soft Handle & Compression Resistant	-----	Systemic vascular repair, primarily for femoral popliteal reconstruction.
Straight Thin Wall	-----	Thin With Soft Handle	Systemic vascular repair, primarily for femoral popliteal reconstruction.
Short Taper	-----	High Strength	Arteriovenous conduit for blood access. Stepped taper and short taper configurations minimise the risk of steal syndrome and high cardiac output ⁴ .
Short Taper Central Support	High Strength & Compression Resistant	-----	Arteriovenous conduit for blood access. Stepped taper and short taper configurations minimise the risk of steal syndrome and high cardiac output ⁴ .

Configuration	Graft Characteristic		Typical Applications
	Supported	Unsupported	
Step Taper	-----	High Strength	Arteriovenous conduit for blood access. Stepped taper and short taper configurations minimise the risk of steal syndrome and high cardiac output ⁴ .
Step Taper - Central Support	High Strength & Compression Resistant	-----	Arteriovenous conduit for blood access. Stepped taper and short taper configurations minimise the risk of steal syndrome and high cardiac output ⁴ .

Contraindications

1. Grafts only to be implanted by Vascular Surgeons who are experienced with the specific techniques required by these medical devices.
2. Thin wall grafts are not recommended for axillo-femoral or axillo-bifemoral bypass reconstruction because of the potential for anastomotic disruption during extreme body movement.
3. Unsupported grafts must not be implanted in such areas where total or partial occlusion can occur as a result of the movement of a patient's body.
4. These prostheses should not be implanted in patients who exhibit sensitivity to ePTFE or materials of bovine origin.
5. These prostheses should not be implanted in the Central Circulatory System as defined by Annex IX, 1.7 of the Medical Device Directive 93/42/EEC with the exception of prostheses implanted using the Blalock-Taussig procedure³.

Warnings

1. Grafts with an external support over the entire length are **NOT** recommended for blood access purposes, however it has been demonstrated that ePTFE (Unity) supported grafts maybe used for such applications⁵.
2. For blood access applications, do **NOT** cannulate the portion of the graft with central spiral support, however it has been demonstrated that this is not applicable for ePTFE (Unity) supported grafts⁵.
3. Do **NOT** remove central spiral support. Attempts to remove this support may damage the graft. If damage occurs, discard the graft.
4. The following precautions should be taken: Point (a) is **particularly important** as gelatin sealed ePTFE grafts feature a degree of longitudinal elasticity :-
 (a) *Excessive tension on the prosthesis must be avoided but moderate tension is essential.*
 (b) *When using embolectomy or balloon angioplasty catheters within the lumen of a graft, the inflated balloon size must be carefully*

matched to the inside diameter of the graft. Failure to size the catheter correctly or the over-inflation of the balloon may result in rupture of either the graft or the balloon.

(c) Cut the graft long enough to ensure that no stress exists on the anastomosis. The patient's body mass and likely extremes of their posture must be considered when determining the length of the graft to be implanted, otherwise stress may be placed on the anastomosis. Failure to consider these aspects may cause anastomotic disruption, resulting in excessive bleeding, loss of function or possible amputation of limb and in the worst case death.

(d) The graft must not be filled with blood before it has been passed through the tunnel. Exerting pressure may cause blood components to be passed through the graft wall leading to the possible formation of seroma⁶.

5. Never cut or open a graft to form a surgical patch.

6. The outer ePTFE wrap on wrapped grafts is **NOT** designed to be removed.

Precautions

1. **DO NOT PRECLOT.** No preclotting required.
2. **DO NOT USE BEYOND THE INDICATED EXPIRY DATE.** The gelatin sealant may not meet the design specification after the expiration date because of hydrolytic action.
3. **DO NOT RESTERILISE. FOR SINGLE USE ONLY.** Do not reuse, reprocess or resterilise. Reuse, reprocessing or resterilisation may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in deterioration of health or death of patients. Reuse, reprocessing or resterilisation may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient end-user.

Figure 1

Standard Graft Soaking Procedure- For All Gelatin Sealed ePTFE Grafts

Standard operating room practice should be followed throughout this procedure to ensure that the sterility of the graft, and of any other items involved, is not compromised.

Stage 1

Fill syringe to 200ml with warm (37°C/98°F) sterile saline.

Stage 2

Add saline solution to the container, taking care not to damage the graft. Rock the blister to ensure complete coverage of the graft.

Stage 3

After 1-5 minutes maximum, the graft can be removed- excess solution is allowed to drain before the graft is implanted. The graft shall not be allowed to dry out after soaking.

4. Clamping may damage any vascular prosthesis. Atraumatic clamps, ideally with soft shod jaws, should be used with a minimum application of force on the unsupported section of graft only. Excessive force must be avoided, as it will damage the external support.

5. **DO NOT PUNCTURE THE GRAFTS IN THE SUPPORTED AREA**, however it has been demonstrated that this is not applicable for ePTFE (Unity) supported grafts⁵.

6. Care must be taken not to puncture the opposite side of the graft when suturing.

7. When using end supported and fully supported grafts there will be less gelatin coverage under the external support. After removing the support the graft could behave more like unsealed ePTFE. In such instances where suture holes are placed within the peeled area the standard implantation protocol should be deployed, eg. use of haemostatic agents.

8. Use a non-cutting tapered needle with non-absorbable monofilament polypropylene or PTFE suture stitched a suitable distance from the graft edge. To minimize suture line bleeding pull the suture at an angle of 90° to the graft.

9. Store in a clean, dry area not less than 0°C (32°F) and not more than 35°C (95°F).

10. **PROSTHESES MUST BE IMPLANTED WITHIN ONE MONTH AFTER REMOVAL FROM THE FOIL POUCH.**

11. Complications which may occur in conjunction with the use of any vascular prosthesis include but are not limited to: redundancy; infection; ultrafiltration or perigraft seroma; thrombosis; embolic events; occlusion; stenosis; mechanical disruption or tearing of the suture line, graft and/or host vessel; excessive suture hole bleeding; formation of pseudoaneurysms due to excessive, localized, or large needle punctures; or perigraft haematomas.

Operative Technique

When implanting a Vascutek Ltd prosthesis, select a bullet tipped tunneller suitably sized to ensure that the graft will have a snug fit. Wetting the graft in a sterile saline solution will facilitate the passage of the graft through the tunneller and optimize the handling characteristics. Grafts should be immersed in saline for no longer than 5 minutes and shall not be allowed to dry out after soaking. See Figure 1 for recommended method.

This is not required for prosthesis rinsed in Rifampicin. If rinsing a prosthesis in Rifampicin prior to implantation Vascutek Ltd recommends following the method as described in Figure 3.



Thrombectomy

Should a postoperative occlusion occur, declotting can be performed by carrying out a longitudinal incision with stay sutures for thrombectomy procedures. When a transverse incision is used, a horizontal mattress suture technique may be required to aid closure.

Spiral Support

To remove the external spiral support from end and fully reinforced Vascutek ePTFE grafts, gently lift the end of the support with gloved hands or atraumatic instruments. Unwind the support by slowly peeling at a right angle to the graft (Figure 2a). Rapid unwinding and/or removal of the support parallel to the axis of the graft may damage the product (Figure 2b). In addition, the use of surgical blades or sharp instruments to remove the support may damage the graft.

Note: Centrally supported grafts have approximately 6cm of spiral reinforcement. This has **NOT** been designed for removal.

to use. Immediate use may increase the risk of haematoma formation.

- The blood access needle should be inserted at an angle of 45° (bevel upward) until the graft has been penetrated. Care must be taken not to puncture the opposite side of the graft.
- Do not cannulate at the same site. Rotate areas to reduce chance of graft damage and the possibility of perigraft haematoma or pseudoaneurysm formation. If using a support graft do not cannulate near or within the supported portion, however it has been demonstrated that this is not applicable for ePTFE (Unity) supported grafts⁵.
- Do not cannulate within one inch of proximal and distal anastomosis.
- To minimize infection, strict adherence to aseptic techniques is recommended.
- Once the needle has been withdrawn, apply moderate digital compression to the cannulation site to assist haemostasis.

Sterilisation

Vascutek Ltd gelatin sealed ePTFE grafts are sterilised in ethylene oxide, supplied sterile, and must not be re-sterilised. The Tyvek® seal on both intermediate and inner trays must be intact. Any damage to the trays renders the prosthesis non-sterile. In the event of damage to the primary packaging, the product must not be used and should be returned immediately to the supplier.




Angioaccess

- When implanting an unsupported prosthesis care must be taken to ensure that a subcutaneous tunnel is created of similar size to the graft, such that a snug fit is achieved.
- Best results are achieved by leaving the graft in place for approximately 2 weeks prior

NOTE: The rifampicin/heparin bonding procedure has not been approved by Health Canada.

Figure 3 Rifampicin Bonding (Optional Procedure)

All Vascutek Ltd gelatin sealed ePTFE grafts may be bonded with the antibiotic Rifampicin using the simple 5 stage procedure outlined below:

Standard operating room practice should be followed throughout this procedure to ensure that the sterility of the graft, and of any other items involved, is not compromised.

Packaging

Trays are enclosed in a foil pouch that serves as a vapour barrier and preserves optimal prosthesis characteristics. A sachet containing a dessicant is included to further aid this purpose.

Note: The foil pouch and outer tray are not sterile. **Only the innermost tray may be introduced to the sterile field.**

Additional Labels

Additional labels are attached for inclusion in patient records to enable the tracking of this device.

Additional Information

Do not expose ePTFE products to temperature greater than 500°F (260 °C). PTFE decomposes at elevated temperatures, producing toxic decomposition products

Disclaimer of Warranties

Many factors are outside Vascutek Ltd's supervision and control after sale of this device. Vascutek Ltd has no control over the conditions under which the device is used, the diagnosis of the patient or the methods or procedures used for implantation. Therefore, Vascutek Ltd makes no warranty or guarantee of this device, expressed or implied.

Vascutek Ltd makes no warranty of merchantability or fitness for a particular purpose. Any warranty or representation by any other person or firm is void. Vascutek Ltd neither assumes, nor authorizes any other person to assume for it any other liability in connection with sale of this device. Vascutek Ltd will not be liable for any incidental or consequential loss, damage, or expense arising directly or indirectly from the use of this device.

Some jurisdictions do not allow the exclusion of or limitations on an implied warranty. Similarly, some jurisdictions do not allow the exclusion or limitations of incidental or consequential damages. Therefore, some of the above exclusions may not apply. This warranty gives specific legal rights. The patient may also have other rights which vary from jurisdiction to jurisdiction.

References

1. Gelatin hydrolysis data on file at Vascutek Ltd.

2. Cavallaro A et al. Pretreatment with gelatin of Dacron Grafts. Experimental Research. Policlinico Sez Chir 1986; 93: 275-283.

3. Tekdogan, M. (1964) Modified Blalock shunt operation in tetralogy of Fallot. Turk J Pediatrics 6:3 129-35.

4. Rosental JJ et al. Prevention of high flow problems of arteriovenous grafts. Development of a new tapered graft. Am J Surg. 1980 Aug;140(2):231-3.

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1. Gelatin hydrolysis data on file at Vascutek Ltd.
2. Cavallaro A et al. Pretreatment with gelatin of Dacron Grafts. Experimental Research. Policlinico Sez Chir 1986; 93: 275-283.
3. Tekdogan, M. (1964) Modified Blalock shunt operation in tetralogy of Fallot. Turk J Pediatrics 6:3 129-35.
4. Rosental JJ et al. Prevention of high flow problems of arteriovenous grafts. Development of a new tapered graft. Am J Surg. 1980 Aug;140(2):231-3.
5. Data on file at Vascutek.
6. Blumenberg RM et al Perigraft seromas complicating arterial grafts Surgery 1985, 97:2 194-204.

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Vascutek Limited: